digital ear thermometer

INSTRUCTIONS FOR USE

(please read and keep for future reference)


tommeetippee.com

EN

digital ear thermometer

INSTRUCTIONS FOR USE

WARNINGS AND PRECAUTIONS

- Correct installation of the hygiene cover ensures accurate readings.
- **Using the thermometer without a hygiene cover may cause inaccurate readings.**
- **Do not submerge the device in any liquids or expose it to direct moisture.**
- This thermometer is not intended as a substitute for a doctor's consultation.
- If the condition does not improve consult your doctor.
- This is not an AP or APG product.
- Intended use: The infrared ear thermometer can measure the infrared heat generated by eardrum area and surrounding tissues to reflect patient's body temperature accurately.
- **Choking from swallowing hygiene covers and connection ring by children or pets is possible, please keep out of reach of children and pets.**

FUNCTIONS

TEMPERATURE ALARM

- If the thermometer detects a body temperature 99.5°F (37.5°C), one long beep will sound followed by three short beeps to warn the user of potential fever.

SWITCHING BETWEEN FAHRENHEIT (°F) AND CELSIUS (°C)

- In 'Power off' mode, press and hold the 'Scan' button.
- Then press the 'On/Mem' button for 3 seconds whilst still holding the scan button. The '°C' will change to '°F'. You can use the same process to change from '°F' back to '°C'.


MEMORY FUNCTION

- Press the 'On/Mem' button again to see the last temperature stored. The thermometer will automatically save the last reading into the memory, if it is within the temperature range of 93.2°F to 108.0°F (34°C to 42.2°C).


INSTALLING THE HYGIENE COVER

Before first use, pull to remove the tab from the battery cover.


Correct installation of the hygiene cover ensures accurate readings. Always use a new and undamaged hygiene cover before each reading.


Gently twist the probe cap to remove.

Place a new hygiene cover into the connection ring.

Insert the thermometer probe into the centre of the hygiene cover and connection ring.


Push until it is secure. If the hygiene cover did not connect firmly, this icon  will be flashing on the LCD screen. Remove the hygiene cover and re-install.


TAKING A READING

HANDY TIP

Ensure the ear canal is clean.


Press 'On/Mem' button. The thermometer is ready for use when the ear icon  stops flashing and two short beeps have sounded.

Gently pull the baby's ear back, straightening the ear canal. Position the probe snugly into the ear canal, ensuring an accurate reading.

To measure the ear temperature, press the 'Scan' button until there is a long beep.

The LCD screen will light up automatically after pressing the scan button, and will automatically turn off after 5 seconds.

The temperature will show clearly on the LCD display.

To power off, hold down the 'On/Mem' button for approximately 5 seconds, until the LCD screen shows 'o'. The thermometer will automatically shut down after 1 minute if it is however not turned off, helping prolong battery life.

TIPS FOR TAKING AN ACCURATE READING

- Always use a new and undamaged hygiene cover before each reading.
- It is recommended that you measure 3 times in the same ear. If the 3 readings are different, select the highest temperature.
- There is no gender and age limitation for using infrared thermometers.
- Holding the thermometer too long may cause a higher ambient temperature reading. This could make the body temperature reading lower than it actually is.

HOW TO INTERPRET YOUR READING

Understanding the meaning of your child's temperature.

READING

>96.3° - 99.4°F (>35.7° - 37.4°C)

>99.4° - 101.3°F (>37.4° - 38.5°C)

>101.3° - 108.0°F (>38.5° - 42.2°C)

MEANING

Acceptable

Elevated

Possible High Fever

REMEMBER if you are concerned about your child's health you should see a doctor. Even in the absence of fever those who exhibit a normal temperature may still need to receive medical attention.

CARE AND CLEANING

The probe is the most delicate part of the thermometer. When cleaning the lens handle with care to avoid damage. If device is accidentally used without a hygiene cover, clean the probe as follows:

1. Use a cotton swab (Alcohol 70% concentration) to clean the lens found inside the probe.
2. Allow the probe to fully dry for at least 1 minute.

- The thermometer should be stored at temperature between -20°C to +50°C.
- Keep the unit dry and away from any liquids and out of direct sunlight.
- The probe should not be submerged in any liquids.
- If the thermometer is used according to the instructions, periodic re-adjustment is not required. If at any time you question the accuracy of the temperature readings please contact us.

TROUBLESHOOTING

Device stabilization in process.

Solution: Wait until  stops flashing.

Battery is low and no more readings are possible.

Solution: Replace the battery.

Reading before device stabilization.

Solution: Wait until  stops flashing.

The ambient temperature is not within the range between 50°F ~104°F (10°C and 40°C).

Solution: Allow the thermometer to rest in a room for at least 30 minutes at room temperature: 50°F ~104°F (10°C and 40°C).

Error 5~9, the system is not functioning properly.

Solution: Unload the battery, wait for 1 minute and then turn on.

Temperature taken is higher than 108.0°F (42.2°C).

Solution: Check the integrity of the hygiene cover and take a new temperature reading.

Temperature taken is lower than 93.2°F (34°C).

Solution: Make sure the hygiene cover is clean and take a new temperature reading.

Device can not be powered on to the ready stage.

Solution: Insert a new battery.

Clinical repeatability 0.18°C newborn to 1 year old, 0.17°C 1-5years old, 0.15°C 5years+

This product comes with a 12-month warranty.

Manufacture Date as the serial number (please open the battery cover, it is shown on the inside of the device.) Ex.SN:E912A000001, the first "E" is External, the second number "9" is the manufacture year 2009, the third and the fourth number "12" is the manufacture month, the remaining figures are the serial number. Our batch code is XXXXVT.

BF type applied part.

CHANGING THE BATTERY

This device is supplied with one lithium cell battery (CR2032).

1. Unscrew battery cover.
2. Flip the battery out with a small screw driver.
3. Insert the new battery under the metal hook on the left side and press the right side of the battery down firmly
4. Replace the battery cover.

The positive (+) side up and the negative (-) side pointed down.

PRODUCT SPECIFICATION

- ☑ Temperature measurement range: 93.2-108°F (34-42.2°C).
- ☑ Operating temperature range: 50-104°F (10-40°C), RH<85%.
- ☑ Storage temperature range: -20-+50°C, RH 85%.
- ☑ Transportation temperature: < than 70°C, RH 95%.
- ☑ Atmospheric pressure: 800~1013 hPa.
- ☑ Enclosure Rating: IP22.
- ☑ Dimensions: 149.1 x 42.5 x 54.9 mm.
- ☑ Weight: 79.7 grams including battery.
- ☑ Battery life: around 3,000 readings. Enable normal use: 1 year.
- ☑ Expected Service Life: 4 years.
- ☑ Complies with ASTM E1965-98, EN12470-5:2003 Clinical thermometers-Part 5: Performance of infra-red ear thermometers (with maximum device), IEC/EN60601-1-2(EMC),IEC/EN60601-1(Safety) standards.
- ☑ Accuracy: 0.4°F (+/-0.2°C) in a temperature range of 95.9-107.6°F (35.5-42°C). Accuracy of +/-0.3°C (0.5°F) for other temperature ranges.
- ☑ This thermometer converts the ear temperature to display its "oral equivalent." (according to the result of the clinical evaluation).

SYMBOL DESCRIPTIONS

Do not reuse.

Please read the instructions for use.

CAUTION.

Manufacturer.

BF type applied part.

IP22

Classification for water ingress and particulate matter.

MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

The digital ear thermometer is intended for use in the electromagnetic environment (for home healthcare) specified below. The customer or the user of the digital ear thermometer should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment-guidance (for home healthcare environment)
RF emissions CISPR 11	Group 1	The digital ear thermometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The digital ear thermometer is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

MANUFACTURER'S DECLARATION-ELECTROMAGNETIC IMMUNITY

The Digital Ear Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Digital Ear Thermometer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Radiated RF IEC 61000-4-3	10 V/m 80 MHz –2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz –2,7 GHz 80% AM at 1 kHz	Recommended separation distance $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range b. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Digital Ear Thermometer is used exceeds the applicable RF compliance level above, the Digital Ear Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Digital Ear Thermometer.

MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The digital ear thermometer is intended for use in the electromagnetic environment (for home healthcare) specified below. The customer or the user of the digital ear thermometer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance (for home healthcare environment)
Electrostatic discharge (ESD) IEC 61000-4-2	Contact±8 kV Air±2 kV, 4 kV,±8 kV,±15 kV	Contact±8 kV Air±2 kV,4 kV,±8 kV,±15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	The digital ear thermometer power frequency magnetic fields should be at levels characteristic of a typical location in a typical home healthcare environment.

RECOMMENDED SEPARATION DISTANCE BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE DIGITAL EAR THERMOMETER

The digital ear thermometer is intended for use in an electromagnetic environment (for home healthcare) in which radiated RF disturbances are controlled. The customer or the user of the digital ear thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the digital ear thermometer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m	
	80 MHz to 800 MHz	800 MHz to 2,7 GHz
	$d = 1,2 \sqrt{P}$	$d = 2,3 \sqrt{P}$
0,01	0,12	0,23
0,1	0,38	0,73
1	1,2	2,3
10	3,8	7,3
100	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

MANUFACTURER'S DECLARATION-ELECTROMAGNETIC IMMUNITY TEST SPECIFICATIONS FOR ENCLOSURE PORT IMMUNITY TO RF WIRELESS COMMUNICATIONS EQUIPMENT

The digital ear thermometer is intended for use in the electromagnetic environment (for home healthcare) specified below. The customer or the user of the digital ear thermometer should assure that it is used in such an environment.

Test frequency (MHz)	Band a)(MHz)	Service a)	Modulation b)	Maximum power (W)	Distance (m)	MMUNITY TEST LEVEL (V/m)	Compliance LEVEL (V/m) (for home healthcare)
385	380 - 390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27	27
450	430 - 470	GMRS 460,FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28	28
710 745 780	704 - 787	LTE Band 13, 17	Pulse modulation b) 217 Hz	0,2	0,3	9	9
810 870 930	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0,3	28	28
1 720 1 845 1 970	1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0,3	28	28
2 450	2 400 - 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0,3	28	28
5 240 5 500 5 785	5 100 - 5 800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0,2	0,3	9	9

NOTE 1 If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.


b) The carrier shall be modulated using a 50% duty cycle square wave signal.

c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

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Copyright © Mayborn USA Inc, 2019
Mayborn USA Inc. dba Mayborn Group