
100 mm

Instructions for use

inearTM

infrared thermometer

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155 mm

Here for you

Your baby may get a fever from time to time and keeping track of their temperature is important.

Our digital in-ear thermometer is quick, accurate and easy to use.

And that's when you can count on us.

Questions?

For further information or support, please scan below to visit our consumer support page at tommeetippee.com

Or join the conversation

EN | IMPORTANT WARNINGS


These instructions can also be found at tommeetippee.com

- This thermometer is not intended as a substitute for consultation with your doctor.
- Adult use only
- Keep out of reach of children when not in use
- Do not allow children to take their own temperature unattended.
- To avoid risk of choking from swallowing small parts or batteries please keep the device out of the reach of children and pets.
- The device should not be submerged into water or other liquids (not waterproof).
- The probe lens is the most delicate part of the thermometer. Use with care when cleaning the probe lens to avoid damage. See care and cleaning section.
- Do not use the thermometer if there are signs of damage to the lens or the thermometer itself. If damaged do not attempt to repair the product.
- No modification of this equipment is permitted.
- Intended for household use only.
- Do not use in direct sunlight
- If the thermometer is used according to the instructions, periodic re-adjustment is not required.

1. Cleaning and storage

The lens in the probe is the most delicate part of the thermometer. When cleaning the lens, handle with care to avoid damage.

- Make sure the scanner lens is clean to ensure an accurate reading.
- Use a 70% alcohol swab or cotton wool moistened with 70% alcohol to gently wipe the lens clean.
- Allow the lens to fully dry for at least 1 minute.
- Never insert a sharp object into the scanner area or any other open surface on the thermometer, because this will cause damage and affect functionality.
- Keep the thermometer dry and away from any liquids and direct sunlight.
- The thermometer should be stored at temperature between -20°C to +50°C.

1. Probe cover
2. Light indicator
3. LCD display
4. Power button / Measure button
5. Probe (take off the probe cover when measuring temperature)
6. Memory / Mute button
7. Mode button (Object temperature, ear temperature)
8. Battery cover
9. Unit switch button (°C/°F)

2. Functions

1. Object temperature mode
2. Ear temperature mode
3. Mute / un-mute
4. Temperature (°C/°F)
5. Low battery
6. Memory recall
7. Temperature value

3. Before first use


Remove the paper slip from the battery compartment, by pulling in the direction shown.


4. How to use - Ear Temperature Mode

A. Remove the probe cover

B. Switch the thermometer on by pressing the power button

C. The thermometer will automatically appear in ear mode

D. Insert the temperature probe into a proper position in the ear canal. Gently pull baby's ear upwards to straighten the ear canal for an accurate reading.

E. Press the measure button to take the temperature

F. The reading will appear on the screen

5. Tips for taking an accurate reading

- Child must be inside for 30 minutes before taking a measurement.
- The child and device must be in the same stable ambient (room) temperature for 30 minutes before operating.
- Always make sure the scanner lens is clean and undamaged.
- Child should not drink, eat or be physically active before / while taking the measurement.
- Remove hats and wait 10 minutes before taking a measurement.
- Don't take a measurement during or immediately after nursing a baby.
- Wait at least 60 minutes after bathing to take a reading.

6. How to interpret your reading

Child (from 3 months+)

Green – OK – 34.0°C – 37.5°C

Red – Fever – 37.6°C – 43°C


7. Memory function

The memory function is great if you want to look back on your last 20 readings, maybe to show a healthcare professional.

When the product is switched on, press the memory button. The screen will display "01" followed by the recorded reading. Press the memory button again to retrieve the second reading, and so on for the past 20 readings

8. How to use – Object Temperature Mode

- A. Remove the probe cover
- B. Switch the thermometer on by pressing the power button
- C. Switch to the object measurement mode by pressing the Mode button 

- D. Point the Infrared sensor to the centre of the object 1-3 cm away, then press and release the Measure button. The object temperature will be displayed on the screen


9. Troubleshooting

	The ambient temperature is higher than 40°C (104°F) or lower than 10°C (50°F). Gently pull baby's ear upwards to straighten the ear canal for an accurate reading.	A long beep, the indicator is red.
	An error occurs when the data is being read from or written to the memory, or the temperature correction is not complete. Press the measure button again and new measurement should be displayed on screen.	A long beep, the indicator is red.
	When the battery voltage is lower than 2.5V ± 0.1V, the low battery symbol will appear on the display. Please replace the batteries.	Silent.

- How to mute – Whilst the product is on, press and hold the memory/mute button for approximately 3 seconds. will appear on LCD display.
- How to resume sound – Whilst the product is on, and the is displayed on the LCD display press and hold the memory/mute button for approximately 3 seconds. A beep will sound and the will disappear.
- Switching between °C/°F: Remove the battery cover, and toggle the °C/°F button in the battery compartment.

10. Replacing the batteries

1. Slide the battery cover off along the marked direction and take it off.
2. Insert the two alkaline AAA batteries into the compartment according to the stated polarities.

- Make sure that the batteries are installed correctly.
- If the low-battery symbol is displayed on the screen, replace the batteries.
- Batteries of the same type should be used. Dispose of used batteries in accordance with the local environmental policies.
- The thermometer is supplied with batteries.
- Not suitable for rechargeable batteries.

FIND OUT HOW TO USE

SCAN QR CODE ON THE BACK OF THE PRODUCT

Product Specifications

Product Name: Infrared Thermometer
Product Model: JPD-FR302
Power Supply Mode: Internal power supply
Operating Voltage: DC 3V
Battery Model: AAA x 2
Operating Mode: Continuous operating
Display: Segment LCD
Measure time: About 1 second
Latency Time: About 3 seconds
Measuring Range: Ear mode: 34.0°C–43.0°C (93.2°F–109.4°F)
Object mode: 0.0°C–100.0°C (32.0°F–212.0°F)
Accuracy (Laboratory): Ear mode: $\pm 0.2^{\circ}\text{C}$ (36.0°C – 39.0°C); $\pm 0.3^{\circ}\text{C}$ (34.0°C – 36.0°C / 39.0°C – 43.0°C); Object mode: $\pm 1.0^{\circ}\text{C}$ / $\pm 2.0^{\circ}\text{F}$
Accuracy (Clinical): $\pm 0.3^{\circ}\text{C}$ ($\pm 0.6^{\circ}\text{F}$)
Measuring site: Ear canal
Reference body site: Armpit
Resolution: 0.1°C (0.1°F)
Mode of operation: Adjusted mode
Memory: 20 temperature readings
Low-battery Alert: The low-battery symbol is displayed if the power voltage is lower than 2.5 V $\pm 0.1\text{V}$
Automatic Power-off: The thermometer automatically powers off if it is not used in 10 ± 1 seconds.
Outer dimensions (mm): 155.5 $\times$ 37.3 $\times$ 39.1 mm
Weight (g): Thermometer (with batteries): 94g
Manufacturing date: see the label
Service life: 2 years
Battery Life: Alkaline dry battery for around 20,000 measurements

Operating Environment:

Temperature: 10°C – 40°C (50°F – 104°F) Humidity: 15%–95% RH, non-condensing Atmospheric pressure: 86–106 kPa
The infrared thermometer has been tested and conforms to the standard ASTM E1965-98. ASTM laboratory accuracy requirements in the display range of 96.8°F to 102.2°F (36°C – 39°C) for ear canal IR thermometers is $\pm 0.4^{\circ}\text{F}$ ($\pm 0.2^{\circ}\text{C}$). Note that for mercury-in-glass and electronic thermometers, the requirement per ASTM Standards E667-86 and E1112-86 is $\pm 0.2^{\circ}\text{F}$ ($\pm 0.1^{\circ}\text{C}$)

Security Class

Type of protection against electric shock: internally powered equipment. Degree of protection against electric shock: Type BF applied part. • Degree of protection against ingress of water: IP22 • Safety degree of using in flammable anesthetic gas blending with air, oxygen or nitrous oxide: Non-AP/APG • No application parts of the thermometer prevents defibrillation charge effect. • No application parts of the thermometer output signal. • The thermometer is permanent installed device.

Storage and Transportation:

The thermometer can be transported using general transportation tools. Severe vibration, shock, or rain must be avoided during transportation. The thermometer must be packaged and then stored in a well-ventilated room without corrosive gas. The ambient temperature must be between -20°C and $+50^{\circ}\text{C}$ (-4°F – 122°F), the relative humidity must be 15%–95% R.H. (non-condensing), and the atmospheric pressure must be 50–106 kPa.

EMC Information-Guidance and Manufacture's Declaration

1* WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

2* WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

3* WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ME equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Declaration - electromagnetic emission	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable

Declaration - electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 8\text{ kV}$ contact $\pm 2\text{ kV}$, $\pm 4\text{ kV}$, $\pm 8\text{ kV}$, $\pm 15\text{ kV}$ air	$\pm 8\text{ kV}$ contact $\pm 2\text{ kV}$, $\pm 4\text{ kV}$, $\pm 8\text{ kV}$, $\pm 15\text{ kV}$ air
Electrical fast transient/burst IEC 61000-4-4	$\pm 2\text{ kV}$ for power supply lines $\pm 1\text{ kV}$ for input/output lines	Not applicable
Surge IEC 61000-4-5	$\pm 0.5\text{ kV}$, $\pm 1\text{ kV}$ line(s) to lines $\pm 0.5\text{ kV}$, $\pm 1\text{ kV}$, $\pm 2\text{ kV}$ line(s) to earth	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Declaration - electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz	Not applicable
Radiated RF IEC 61000-4-3	10V/m 80 MHz to 2.7 GHz	10V/m

Declaration - IMMUNITY to proximity fields from RF wireless communications equipment					
Immunity test	IEC 60601 test level				Compliance level
	Test frequency	Modulation	Maximum power	Immunity level	
Radiated RF IEC 61000-4-3	385 MHz	**Pulse Modulation: 18Hz	1.8W	27 V/m	27 V/m
	450 MHz	*FM+ 5Hz deviation: 1kHz sine	2 W	28 V/m	28 V/m
	710 MHz 745 MHz 780 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m
	810 MHz 870 MHz 930 MHz	**Pulse Modulation: 18Hz	2 W	28 V/m	28 V/m
	1720 MHz 1845 MHz 1970 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m
	2450 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m
	5240 MHz 5500 MHz 5785 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m

Note* - As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

Note** - The carrier shall be modulated using a 50 % duty cycle square wave signal.

CAUTION. **MD** Medical Device  BF type applied part.

IP22

Classification for water ingress and particulate matter.

Indicates this device is subject to the Waste Electrical and Electronic Equipment Directive in the European Union. To protect the environment dispose of useless device at appropriate collection sites according to national or local regulations.

Please read instructions.  Battery recycling.

Information on the Disposal of Waste Electrical and Electronic Equipment (WEEE)

This symbol means that your product shall be disposed of separately from household waste in accordance with local laws and regulations. When this product is ready for disposal, take it to a collection point designated by local authorities. Separating and recycling of your product at the time of disposal will help conserve natural resources and ensure that the product is recycled in a manner that protects human health and the environment. Penalties may be applicable for the incorrect disposal of this product in accordance with national legislation.

Manufactured by: Shenzhen Jumper Medical Equipment Co., Ltd.
Address: D Building, No. 71, Xintian Road, Fuyong Street,
Baoan, Shenzhen, Guangdong, China. 518103.

Distributed by: Mayborn (UK) Limited, Balliol Business Park,
Newcastle upon Tyne, NE12 8EW.
Mayborn France Sarl, 92100, France.
Mayborn USA Inc. CT 06901, USA.

SUNGO Certification Company Limited
3rd floor, 70 Gracechurch Street, London. EC3V 0HR.

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生產廠商: Shenzhen Jumper Medical Equipment Co., Ltd. Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China, 518103.

Distributed by: Mayborn (UK) Limited, Balliol Business Park, Newcastle upon Tyne, NE12 8EW.
Mayborn France Sarl, 92100, France.
Mayborn USA Inc. CT 06901, USA.

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Shenzhen Jumper Medical Equipment Co., Ltd. Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China, 518103.

Distributed by: Maybom (UK) Limited, Balliol Business Park, Newcastle upon Tyne, NE12 8EW.
Maybom France Sarl, 92100, France.
Maybom USA Inc. CT 06901, USA.

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TOMMEE TIPPEE®. IN EAR™. Copyright © 2024.
Mayborn (UK) Ltd
All IP rights held by Mayborn (UK) Limited