
**tommeetippee**<sup>®</sup>

# nasal aspirator

instructions for use

model number: NC003

[tommeetippee.com](http://tommeetippee.com)

0423773

## **EN | instructions for use**

Please read and retain for future reference. These instructions can also be found at **[tommeetippee.com](https://tommeetippee.com)**

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
## EN | PARTS DESCRIPTION

- ① Removable cover for disposable cups. Twist to remove
- ② Silicone ring
- ③ On / Off button
- ④ Battery cover
- ⑤ Disposable cups x 25
- ⑥ Compartment for disposable cups
- ⑦ Pink silicone tip 0-3 Months
- ⑧ Blue silicone tip 3-12 Months
- ⑨ Orange silicone tip 12 Months +
- ⑩ Spare silicone sealing ring


Any questions?  
Visit our website: [tommeetippee.com](http://tommeetippee.com)

## INSERTING / CHANGING BATTERY

1. Slide off battery cover.
2. Insert 2 x AA alkaline batteries.
3. Slide battery cover back into position until you hear a click.  
Do not use rechargeable batteries.

## HOW TO USE

1. Clean hands before use.
2. Wash and Sterilise clear cover, sealing ring and aspirator tip. Disposable cup can be re-used if handwashed thoroughly between use.
3. Place disposable cup inside unit, position sealing ring in place then position clear cover over the top – twist to seal and lock in place.
4. Select the correct tip for age of baby / toddler.
5. Ensure baby's mouth is open.  
Place tip into baby nose. DO NOT place deeply into the nose.  
DO NOT press the tip hard against the nasal wall.  
DO NOT apply the device to the ear, eye or throat this product is designed to remove nasal mucus only.  
When using the aspirator DO NOT block baby's other nostril.
6. Press "i" to operate the unit only in bursts of a few seconds.
7. Use for a few seconds only to remove mucus.
8. Turn off the device by pressing "O".
9. Disassemble cover and remove disposable cup holding the nasal mucus.  
Dispose of mucus.
10. Wash, sterilise, dry the cover and silicone tip. Wash and dry the disposable cup or use replacement provided. Reassemble cover, disposable cup and silicone tip.
11. Wipe the main body with soft cloth. DO NOT immerse the unit in water.  
Always use unit with disposable cup to collect mucus.  
Always change or empty disposable cup of mucus when the cup becomes two thirds full.  
Do not turn aspirator upside down when mucus is in the cup as it may block the aspirator.

## CARE AND USE

When using the aspirator make sure the baby opens their mouth.

When using the aspirator do not block baby's other nostril.

Clean the aspirator after each use.

Do not store the unit with mucus inside the disposable cup. Always clean / sterilise after each use.

Please remove batteries if the aspirator will be left unused for longer than 14 days.

Always insert silicone sealing ring correctly before each use.

Keep out of reach of children.

Do not use the aspirator if there is any wound, infection or bleeding of the nose and consult your doctor.

Only use the accessories provided with the aspirator.

## TROUBLE SHOOTING

Problem	Solution
No aspiration	Check silicone sealing ring is in position Check cover is locked and sealed into place Check Tip is correctly assembled
On / Off button not operating	Check batteries are inserted correctly Change batteries

## GUIDELINES

This device complies with the EU Directive 93/42/EEC concerning medical products, and the European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility.

## CLASSIFICATION

1. Internally powered equipment;
2. Type BF applied part;
3. IP22;
4. Sterilised or disinfection: Steam steriliation;
5. Category AP / APG equipment: N/A.
6. Continuous operation.

## EXPLANATION OF SYMBOLS

Symbol	Description of symbol
	Application part, Type BF
	Direct current
IP22	Protected against access to hazardous parts with a finger and against vertically falling water drops when enclosure tilted up to 15°
	Refer to instruction manual / booklet
	DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary
	Manufacturer
	Serial number
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical device
	Authorized representative in the European Community/European Union
	Importer

## EMC DECLARATION

The ME EQUIPMENT (NASAL ASPIRATOR) or ME SYSTEM is suitable for home healthcare and clinic environments.

**WARNING:** Do Not use this equipment near active HF surgical equipment or the RF shielded room of an EM magnetic resonance imaging system, where the intensity of EM disturbances is high.

**WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

**WARNING:** Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

**WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm (12 inches) to any part of the equipment NCO03, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

1. All necessary instructions of maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE regarding electromagnetic disturbances for the expected service life.  
Portable and mobile RF communication equipment may affect the performance of the NCO03, avoid strong electromagnetic interference when in use, such as near mobile phones, microwave ovens, etc.
2. Guidance and manufacturer's declaration - electromagnetic emissions and immunity.

## GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS AND IMMUNITY

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable

Table 2

<b>Guidance and manufacturer's declaration - electromagnetic immunity</b>		
<b>Immunity test</b>	<b>IEC 60601-1-2 Test level</b>	<b>Compliance level</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV ±4 kV, ±8 kV ±15 kV air	±8 kV contact ±2 kV ±4 kV, ±8 kV ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC 61000-4-6	Not applicable	Not applicable
Radiated RF IEC 61000-4-3	10 V/m 80MHz - 2,7 GHz 80 % AM at 1 kHz	10 V/m 80MHz - 2,7 GHz 80 % AM at 1 kHz
NOTE UT is the a.c.mains voltage prior to application of the test level.		


Table 3

Guidance and manufacturer's declaration - electromagnetic immunity							
Radiated RF IEC 61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
	450	430-470	GMRS 460 FRS 460	FM $\pm$ 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
	745						
	780						
	810	800-960	GSM 800/900 TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
	870						
	930						
	1720	1170-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
	5500						
	5785						

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TOMMEE TIPPEE®. Copyright © 2021.
**Guangzhou Berrcom Medical Device Co., Ltd.**
**Shanghai International Holding Corp. GmbH (Europe)**
**Mayborn (UK) Limited**, Balliol Business Park,
(UK) Limited and/or its affiliates.